

**scan**



**CLARiTI**

ADRC Consortium for Clarity in ADRD Research Through Imaging

# **PET Technical Procedures Manual**

# PET Technical Procedures Manual

## Document Version History

Version	Date	Author	Description
1.12	Jan 21, 2022	Suzanne Baker	Replaced QuickStart Guide with PET checklist
1.2	July 22, 2022	Suzanne Baker	Added GE Disc ST reconstruction sheet
1.21	March 9, 2023	JiaQie Lee, Suzanne Baker, Brittany Hale	Adopted consistent language for Subject ID, PTID & NACCID; Updated PET Upload Form
1.22	April 5, 2023	Trevor Chadwick	Revised FDG and PI2620 acquisition times
1.23	June 22, 2023	Trevor Chadwick	Revised SCAN PET Checklist; Step 5 of Appendix D
1.24	May 6, 2024	Trevor Chadwick	Added GE Omni Legend PET/CT scanner, revised GE Discovery MI and GE Sigma scanner appendices
1.25	June 17, 2024	Trevor Chadwick	Added radiotracer Flutemetamol and preferred acquisition time for NAV4694, PI2620, and FDG
1.26	August 15, 2024	Wesley Thomas, Suzanne Baker	Updated Discovery 690 XT parameters Removed MK6240 and PI2620 non-preferred acquisition times. Added "allowable" for PIB and FDG non-preferred acquisition times. Included Clariti logo.
1.27	April 25, 2025	Trevor Chadwick	Added radiotracer RO948

# PET Technical Procedures Manual

## TABLE OF CONTENTS

<b>SCAN PET Checklist</b>	<b>5</b>
<b>General Information</b>	<b>7</b>
<b>Contact Information</b>	<b>7</b>
<b>Site Qualification</b>	<b>7</b>
Hoffman Phantom Scans .....	7
<b>Regulatory</b>	<b>9</b>
<b>Quality Monitoring of PET scanner</b>	<b>9</b>
PET-only Scanner .....	9
PET/CT Scanner .....	9
PET/MR Scanners .....	9
<b>PET General Information</b>	<b>10</b>
Ambient Conditions .....	10
Naming convention .....	10
Documentation .....	11
<b>PET Procedures</b>	<b>12</b>
Radiotracer Doses .....	13
Preparing subject .....	14
Drawing up, assaying, and injection dose .....	14
Subject Positioning .....	15
Attenuation Correction .....	16
Emission Acquisition .....	16
Post PET Scan .....	17
<b>Appendix A: PET Scanner Reconstruction Parameters</b>	<b>18</b>
GE Discovery MI - PET/CT scanners .....	19
GE Discovery 600, 610, 690, 710, and MI DR - PET/CT scanners .....	20
GE Discovery 690-XT - PET/CT scanner .....	21
GE Discovery RX - PET/CT scanners .....	22
GE Discovery STE - PET/CT scanners .....	23
GE Discovery ST - PET/CT scanners .....	24
GE Omni Legend PET/CT scanner .....	25
GE Signa - PET/MR scanners .....	26
Philips Vereos - PET/CT scanners .....	27
Philips Ingenuity TF - PET/CT scanners .....	28
Philips Gemini TF - PET/CT scanners .....	29
Siemens BioGraph Horizon - mCT TrueV PET/CT scanners .....	30
Siemens BioGraph mMR .....	31
Siemens BioGraph Vision PET/CT scanners .....	32
Siemens BioGraph mCT - mCT TrueV PET/CT scanners .....	33
Siemens BioGraph HiRes – 81 slice PET/CT (Model 1080) .....	34
Siemens BioGraph TruePoint – and TruePoint TrueV PET/CT scanners (Models 1093, 1094) ...	35
Siemens ECAT Exact HR+ (BGO) 63-slice scanners .....	36
<b>Appendix B: IDA Uploader Introduction and User Registration (SCAN only)</b>	<b>37</b>
<b>Appendix C: Obtaining/Installing IDA-uploader (SCAN only)</b>	<b>40</b>
<b>Appendix D: How to Upload PET Data to SCAN using IDA-uploader (SCAN only)</b>	<b>49</b>

# PET Technical Procedures Manual

<i>Appendix E: How to Upload PET Hoffman Phantom Data to SCANQUAL (SCAN only)</i>	<b>59</b>
<i>Appendix F: SCAN PET Upload Form (SCAN only)</i>	<b>64</b>

# PET Technical Procedures Manual

**This Technical manual was originally developed for the SCAN project.** The CLARiTI project uses the same PET imaging techniques as the CLARiTI project. Therefore, information about PET scanners, tracers, acquisition times and reconstruction are identical for the 2 studies. Please follow these methods for SCAN and CLARiTI and refer any questions to Suzanne Baker.

## SCAN PET Checklist

About this checklist: This checklist is designed to serve as a roadmap

### Step1. Identify SCAN Liaison(s) Responsible for Study Oversight and Submit Site Information

- Description of SCAN Liaison role:
  - Respond to queries from the SCAN team
  - Maintain list of site contacts for imaging site
  - Ensure protocol compliance
  - Submit PET contact information to SCAN via this form

### Step2. Submit SCAN site information

- Submit Site and Contact Information team via [Site Information Form](#) noting:
  - Who will serve in the PET liaison role
  - All site contacts for imaging site (Technologists/Uploaders/PI's)

### Step3. Prepare to participate in the SCAN Study

- Ensure regulatory compliance and institutional approval for data sharing
- Get ready to collect images:
  - Qualify your PET Scanner
  - Know how to order and allocate correct PET tracers
  - Review PET upload form with imaging center staff responsible for the protocol
- Identify participants that are eligible for the study:
  - All participants must have:
    - A NACC ID (NOTE: If the participant does not have an assigned NACC ID, the data upload will not be accepted)
    - Been scanned after Jan 1, 2021
    - PET data collected according to SCAN standards

### Step4. Collect and upload data according to protocol

- Prepare for image collection:
  - Check-in with imaging center to review the following:
    - Protocol for data acquisition will be used
    - Protocol for image reconstruction will be used

# PET Technical Procedures Manual

- Technician acquires images:
  - Follow protocol for image acquisition tracer and dose
  - Ensure PET upload form is filled out during each scan
  - Reconstruct images as per scanner specifications (Appendix A)
- Upload images:
  - Ensure images are uploaded properly (Appendix D)
  - Make sure you have the filled-out PET upload form
  - After images are uploaded you will transfer PET upload form via web interface
- Troubleshooting:
  - Identify and correct any mistakes that are detected in image reconstruction or upload
  - Be aware of appropriate SCAN personnel available for questions/advice

# PET Technical Procedures Manual

## General Information

The purpose of this manual is to further explain the PET imaging component of the SCAN protocol. Standard procedures are needed to ensure consistency of data collection in this study. Contents include:

- Imaging center training and setup process
- Image acquisition guidance (required image acquisition schedule, tracer, administration and recommended techniques)
- Submission of image data to the SCAN Data Portal

This manual contains information for study-site clinical staff involved with the care of study participants during the imaging procedure and those involved with the processing and transfer of PET imaging data.

## Contact Information

### **Questions about Uploading SCAN data:**

[data.coordinator@loni.usc.edu](mailto:data.coordinator@loni.usc.edu)

(For questions/concerns regarding individual subjects contact the study coordinator at your referral site).

### **Technical/QC questions:**

[koeppe@umich.edu](mailto:koeppe@umich.edu) (Robert A. Koeppe) or [slbaker@lbl.gov](mailto:slbaker@lbl.gov) (Suzanne L. Baker).

For questions regarding scanner specific acquisition and reconstruction parameters.

## Site Qualification

It is preferable for sites to use existing qualified TRC-PAD, ADNI, LEADS, DIAN, DIAN-TU, Pointer, or NiAD scanners for PET imaging. If you are using a scanner that has not been qualified for one of these projects by Bob Koeppe, it will need to be qualified before imaging can be performed. Please contact Bob Koeppe ([koeppe@umich.edu](mailto:koeppe@umich.edu)).

If you plan to acquire data for SCAN on a new PET scanner, please contact Bob Koeppe. You will need to perform two Hoffman phantom scans and send images to Bob Koeppe prior to scanning any subjects on the new PET scanner.

## **Hoffman Phantom Scans**

### **Do you need to perform Hoffman phantom scans?**

Phantom scans will be performed during the site qualification process, to set up SCAN PET acquisition and reconstruction protocols and to validate their quality. It will also be required after a major hardware upgrade of your PET system, to detect any potential deviation due to the upgrade. Once the system is up and running again, please check the potential impact of the upgrade on the acquisition parameters and communicate

# PET Technical Procedures Manual

with Bob Koeppe (koeppe@umich.edu) to find out what change is considered acceptable.

**NOTE:** If you are using a new scanner that is not listed in the appendix, then the University of Michigan PET QC Team will provide you with scanner specific instructions for your site's PET scanner.

## Instructions on Hoffman phantom scans:

### 1. Filling Hoffman:

**Radioactivity:** ~0.6 mCi F-18 at time of scan start.

Fill the phantom carefully to remove air bubbles as much as possible. Mix thoroughly.

Tips for successful filling:

- a) Take out ~50 ml from a filled phantom, BEFORE putting in the radioactivity. It works best to use a spinal needle to inject the radioactivity. Be careful not to “stab” the plastic with the tip of the needle as this often creates a “hot” spot in the scan.
- b) Next, fill the phantom with all but ~10 ml, and mix thoroughly.
- c) Then put in all but the last ~1 ml, and mix again.
- d) Finally, fill that last ml or so as best you can to get rid of the remaining air bubble and mix on last time. For mixing, it works well to put down absorbent pads and have two people roll the phantom back and forth for a half a minute to a minute, which allows less direct contact with the phantom, thus reducing personnel exposure. (One can even use the floor and their feet.)

### 2. Scanning of the Hoffman phantom:

The Hoffman phantom should be positioned in the scanner in the same manner as a subject laying on their back. The frontal cortex should be closer to the top of the scanner, occipital cortex to the floor.

**Acquisition:** The acquisition time for the Hoffman phantom scan should be four 5-min frames, for a total of 20 min. **Reconstruction:** The same reconstruction parameters should be used as the human scans.

### 3. Upload Hoffman phantom reconstruction files:

See [Appendix B](#) part D.

### 4. Phantom scan evaluation and result:

University of Michigan QC will examine the phantom data and determine if the correct parameters have been met and assure there are no other underlying problems with the scanning session. Your site will be notified by email if the phantom scans pass or fail, and whether your PET system has been certified for SCAN.



# PET Technical Procedures Manual

## Regulatory

You must have institutional approval to share your data. Sites must be appropriately licensed through appropriate state or federal agencies to receive and use all radiotracers prior to imaging.

Sites must also receive both IRB approval and radiation safety committee (RSC) or similar approval, before scanning any subjects.

## Quality Monitoring of PET scanner

To ensure scanner/ancillary equipment stability and quality, each site is required to perform ongoing quality control procedures.

### **PET-only Scanner**

- The PET scanner should have an up-to-date calibration and normalization on the date of each imaging session.
- A daily QC/blank scan should be done at the beginning of the day the scanning is to be completed. This scan should be visually inspected for abnormalities. If there is a possibility that the abnormality could impact the quality of the PET scan the study should be rescheduled.

### **PET/CT Scanner**

- The PET/CT scanner should have an up to date calibration and normalization on the date of the imaging session.
- A daily QC check should be done at the beginning of the day the scanning is to be completed. This scan should be visually inspected for abnormalities. If there is a possibility that the abnormality could impact the quality of the PET scan the study should be rescheduled.
- Daily CT should be performed as recommended by the specific vendor, but typically should include a "checkup/calibration" procedure and a water phantom scan. The checkup/calibration procedure guarantees optimum image quality by warming up the x-ray tube and should be performed at startup and within 1 hour prior to any scan. The water phantom provides quality measurements of 3 parameters. The parameters are the CRT value of water calculated in Hounsfield units (HU), the pixel noise of images calculated as a standard deviation, and the tube voltages measured directly on the x-ray tubes. These three measurements should be determined for all available Kvp values.

### **PET/MR Scanners**

- The PET/MR scanner should have an up-to-date calibration and normalization on the date of the imaging session.

# PET Technical Procedures Manual

- On the day of a scan, a daily normalization and detector setup should be completed before any subject scanning is performed. This scan should be visually inspected for abnormalities. If there is a possibility that the abnormality could impact the quality of the PET scan the study should be rescheduled.
- Daily MR QC should be performed as recommended by the specific vendor.

## Ancillary Equipment

- Quality control of the dose calibrator should be performed throughout the course of the study. This typically will include daily constancy, quarterly linearity checks and annual accuracy tests

## PET General Information

### Ambient Conditions

Standardization of the environment during the uptake period following **tracer** administration is not essential.

### Naming convention

It is VERY important that each site follow standard file identification so that all scans can be easily identified.

When uploading data to SCAN, you will enter your PTID or Subject ID (local Subject ID) which is a sequence of up to 10 characters. Formats vary by each center that you use for UDS (uniform dataset) identification internally. As a rule, enter the complete ID. However, you may leave off leading zeros if the format is all numeric digits. Also, either upper or lower case may be used for any alphabetic characters.

The Subject ID or PTID is different from the NACCID.

- Subject ID or PTID - this is the ADRC-managed participant ID, or "local ID". Formats vary by center.
- NACC ID - this is the NACC-managed participant ID. It is a string with the prefix 'NACC' followed by 6 digits.

PLEASE NOTE: If the participant does not have an affiliated NACCID assigned to the Subject ID or PTID you entered, the upload will not be accepted.

The tracer should be included in the Series Description:

In the PET scan Series Description, include PET radiotracer used as follows:

Amyloid\_PIB

Tau\_Flortaucipir

# PET Technical Procedures Manual

Amyloid\_Florbetapir  
Amyloid\_Florbetaben  
Amyloid\_NAV  
Amyloid\_Flutemetamol  
FDG

Tau\_MK6240  
Tau\_PI2620  
Tau\_GTP1  
Tau\_RO948

## Documentation

The study coordinator must ensure the PET Technologist has a copy of all PET Upload Form prior to each scan session, and that the PET scan information form matches the type of study being performed (amyloid, tau, or glucose metabolism). Be sure to complete the PET Upload form as the study is being acquired. A process should be established for transferring this form(s) back to the study coordinator.

# PET Technical Procedures Manual

## PET Procedures

Be sure to have the PET scan upload form and fill out accordingly during the procedure.

**ORDERING TRACER:** Radiotracers are generally ordered from radio-pharmacies by the imaging center. Please be certain that your research participants will receive the correct tracer at the correct dose and be aware of the delivery time because the radioactivity will decay.

Considerations for acquisition and reconstruction:

**1. Reconstruction framing:** Any frames within the Acquisition start-stop time listed in [Table 1](#) should be reconstructed as multiple 5-minute frames.

*Examples: for Florbetaben 90-110, this should be 4x5min frames; for MK6240 70-110 this should be 8x5min frames; for GTP1 this should be 6x5min frames.*

**2. Are longer dynamic scans (emission started at time of injection) acceptable?**

As long as the dynamic scan fully includes at least one of the start-stop times listed in [Table 1](#) for that tracer, this is acceptable. You are encouraged to upload the reconstruction of the full dynamic dataset as long as the framing within the tracer-specific acquisition time is 5 minute frames. Examples: PIB 0-50min would not be acceptable, but a PIB scan from 0-60min would be acceptable where the framing from 0-40min should adhere to your site's reconstruction protocols and framing within 40-60min must be 4x5min frames

**3. Are longer static scans acceptable?**

This is now acceptable (also preferred) for acquisition and reconstruction as long as it includes a preferred acquisition time. Example: If you acquire Flortaucipir from 75-115min post-injection (which is acceptable because it includes the preferred 80-100min acquisition time), we would rather you fill out the PET meta data form with the write-in start and stop times of 75 and 115 and upload 8x5min frames, instead of doing a special reconstruction that is only 80-100 min as 4x5min frames.

**4. We are acquiring something different than what is listed for the acquisition times, what should we do?**

- **If the protocol includes the acquisition start and stop time but collects more data, then whenever possible reconstruct only the preferred time:** ex if you are acquiring PIB 0-70min, we ask you upload the data from 50-70min only
- **If your current protocol does not fully encompass the acquisition start and stop times, then this is no longer considered SCAN compliant.**

# PET Technical Procedures Manual

## Radiotracer Doses

Table 1: Radiotracer Doses, Acquisition times post-injection

	Suggested Target Dose ± 10%		Minimum injectable dose		Acquisition start-stop time post-injection (min)
Tracer	mCi	MBq	mCi	MBq	
AMYLOID TRACERS					
PIB	15	555	8	300	Preferred: 50-70 Allowed: 40-70, 40-60
Florbetapir	10	370	6	225	50-70
Florbetaben	8	295	5	185	90-110
NAV4694	8.1	300	6.5	240	50-70
Flutemetamol	5	185	4	150	90-110
TAU TRACERS					
Flortaucipir	10	370	6	225	80-100
MK6240	5	185	4	150	90-110
PI2620	5	185	4	150	45-75
GTP1	7	260	5	185	60-90
RO948	10	370	5	185	70-90
GLUCOSE METABOLISM TRACER					
FDG	5	185	4	150	Preferred: 30-60 Allowed: 30-45

The above are suggested target doses  $\pm$  10%. We will accept data scanned with the above tracers as long as the amount injected does not fall below the minimum injectable dose and scans are acquired using the start-stop times listed under Acquisition time post-injection. If you are acquiring data using an acquisition time not listed, please contact Suzanne Baker ([slbaker@lbl.gov](mailto:slbaker@lbl.gov)).

# PET Technical Procedures Manual

Using aseptic technique and radiation shielding, draw of the radiotracer and assay with a dose calibrator. Record the assay time to the nearest minute. Do not add saline to the dose prior to administration. Adding saline could potentially lead to precipitation out of solution form.

## Preparing subject

- Have the subject use the restroom and empty their bladder.
- Allow them to lie comfortably in a bed or reclining chair in a room. Supply them with blankets/pillows as needed to maximize their comfort.
- Obtain intravenous access using a small angiocath.

## Drawing up, assaying, and injection dose

- Draw up sufficient tracer to achieve the target dose and assay with a dose calibrator.
- ***Record the assay time to the nearest minute.*** Do not q.s. (add saline) to the dose prior to administration. Adding saline could potentially lead to precipitation out of solution form.
- Inspect the radiopharmaceutical dose solution prior to administration and do not use it if it contains particulate matter or is discolored.
- Inject the radiotracer and follow the injection with an intravenous flush of 0.9% sterile sodium chloride. ***Record the injection time to the nearest minute.*** The IV line can be discontinued at this time.
- Re-assay the dose syringe. If the residual activity is 0.1 mCi or greater, record the amount and correct the amount of the injected dose for the residual activity.
- If acquisition does not start with injection, allow the subject to rest comfortably in the room during the incorporation period until about 10-15 min prior to the start time for the particular radiotracer.
- 10-15 minutes before scanning, have the subject use the restroom and empty their bladder.
- For PET/CT scanners, the subject should be placed on the scanning table

enough to obtain a CT scan prior to the emission scan, and still be able to begin the emission acquisition at the appropriate start time. **PET-only systems will acquire a transmission scan following the emission acquisition**, so preparation of the subject can begin a few minutes later than if on a PET/CT system.

## Subject Positioning

Proper subject positioning is a key aspect of the successful completion of the PET exams. It is important to take the time necessary to ensure not only that the subject is properly positioned but also can **comfortably** maintain that position throughout the duration of the scanning session. **Excessive motion and in particular a difference in the subjects' position between the emission scan and the transmission or CT scan used for attenuation correction is the single most common cause of failed studies.**

- Have the subject remove any bulky items from their pockets such as billfolds, keys, etc.
- In addition, they should remove eyeglasses, earrings, and hair clips/combs if present. If possible, they should try and remove hearing aids also.
- Position the subject so that their head and neck are relaxed. It may be necessary to add additional pads beneath the neck to provide sufficient support. Use the lasers to ensure there is little or no rotation in either plane. The head should be approximately positioned such that the PET scanning planes are parallel to the imaginary line between the external canthus of the eye and the external auditory meatus (orbitomeatal plane) and the head is centered in the sagittal plane. More important exact matching of the orbitomeatal plane is making sure the subject is comfortable, which will hopefully translate to less subject motion during the scan.
- Alignment marks should be put on the subject using the laser system, which can then be subsequently used to check alignment and reposition the subject as necessary.
- Use support devices under the back and/or legs to help decrease the strain on these regions. This also will assist in the stabilization of motion in the lower body.
- Once the subject has been positioned foam pads can be placed alongside

# PET Technical Procedures Manual

the  
head

for additional support. Velcro straps and/or tape should also be used to secure the head position. Vacuum beanbags can also be used in this process.

- The subjects should be offered a “panic button” or be reassured that someone is watching or able to hear them at all times.
- Proper positioning of the subject to get the entire head in the field of view is critical to the success of the project.

Checking the subject positioning and readjusting (if possible) the position of the subjects’ head should be done often throughout the study.

## Attenuation Correction

### PET/CT Scanners

- Standard CT acquisition parameters, but low effective mAs (~30 is typical)
- The subject must undergo the CT scan starting about 5 minutes prior to starting the PET scan. Be sure to prepare the subject so that you are ready to press “start” for the PET emission scan at the required time.

### PET-only Scanners

- For dynamic PET scans (emission starts at injection), acquire the transmission scan for attenuation correction for 5-6 minutes before the emission scan. The subject should be repositioned “on their marks” prior to acquiring the emission scan.
- For static PET scans, acquire the transmission scan using rod sources for 5-6 minutes after the acquisition of the standard emission scan. The subject should be repositioned “on their marks” prior to acquiring the transmission scan.
- Segmentation and re-projection routines will be applied for attenuation correction.

## Emission Acquisition

- Every effort should be made to make sure the acquisition time fully encompasses the target start-stop time from [Table 1](#). Scheduling scans very close in time following a clinical study (which could run late) is not an acceptable reason for starting a scan



# PET Technical Procedures Manual

late.

*Example: If you are only acquiring Flortaucipir data from 80-100min, every effort should be made to start the acquisition as close to 80min post-injection as possible.*

- It is crucial that the subject's position is checked several times throughout the PET scan. A good idea is to check the subject's marks using the laser system at the end of each 5 min scan frame. The subject's position should be returned as closely as possible to the original position just at the beginning of the next scan frame.

## Post PET Scan

- Reconstruct images using parameters specific to your site's scanner (See Appendix A). The framing should always be multiple 5min frames within the start-stop time window in [Table 1](#). Otherwise, the same reconstruction parameters should be used for all emission scans. Upon completion of the reconstruction, review all the images to assess for artifacts.
- Archive ALL raw and processed study data including copies of the attenuation correction-related files (transmission, CT, etc), normalization and blank scans. It is necessary to archive and store raw and processed data at the imaging site for the duration of the project and local guidelines (approximately 5 years).
- Upload PET data using IDA-Uploader (Appendix B)

# **Appendix A:**

## **PET Scanner Reconstruction Parameters**

# PET Technical Procedures Manual

GE

## Discovery MI - PET/CT scanners

(Arizona, Univ of Washington, Mayo)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

**Randoms Correction:** Singles: (not 'delays')

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

#### Reconstruction Method:

**VPHD or VPFX** (not VPHDS or VPFXS). Use VPHD or VPFX, which ever your site normally uses

**6 iterations; 16 subsets**

**Grid: 192 x 192 x (53, 71, or 89)** (depending on the Disc-MI model)

**FOV: 256 mm** (results in voxel size of 1.333 mm)

**Slice Thickness: 2.79 mm**

**Smoothing Filter: NONE or 0.0** (for all filter options: loop filter, post-filter and z-axis filter)

All corrections '**On**'

**Note:** Reconstruction offsets should be used to assure the head is fully in the in-plane FoV.

# PET Technical Procedures Manual

GE

Discovery 600, 610, 690, 710, and MI DR - PET/CT scanners

(Wake Forest, Arizona, Oregon, Kansas, Boston)

## Acquisition Parameters:

### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

### Randoms Correction:

Singles: (not 'real-time subtraction' or 'delays')

## Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

### Reconstruction Method:

**VPHD or VPFX** (not VPHDS or VPFXS). Use VPHD or VPFX, which ever your site normally uses

**4 iterations; 24 subsets (or as close to 24 subsets as the software allows)**

**Grid: 192 x 192 x 47**

**FOV: 256 mm** (results in voxel size of 1.333 mm)

**Slice Thickness: 3.27 mm**

**Smoothing Filter: NONE or 0.0** (for all filter options: loop filter, post-filter and z-axis filter)

All corrections '**On**'

**Note:** Reconstruction offsets should be used to assure the head is fully in the in-plane FoV.

# GE **PET Technical Procedures Manual**

## Discovery 690-XT - PET/CT scanner

(Mayo)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

**Randoms Correction:** Singles: (not real-time subtraction)

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

#### Reconstruction Method:

**VPHD**

**4 iterations; 20 subsets**

**Grid: 256 x 256 x 79**

**FOV: 300 mm** (results in voxel size of 1.17 mm)

**Slice Thickness: 1.96 mm**

**Smoothing Filter: NONE or 0.0** (for all filter options: loop filter, post-filter and z-axis filter)

All corrections '**On**'

**Note:** Reconstruction offsets should be used to assure the head is fully in the in-plane FoV.

This scanner says Discovery 690 in the header, but the 79 (instead of 47) planes and 1.96mm slice thickness tells you it's the XT.

# PET Technical Procedures Manual

GE

Discovery RX - PET/CT scanners

(Hopkins)

## Acquisition Parameters:

### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

### Randoms Correction:

Singles: (not real-time subtraction)

## Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Reconstruction Method: 3D Iterative** (fully 3D Iter; not 3D FORE IR):  
**4 iterations; 21 subsets**

**Grid: 128 x 128**

**FOV: 256 mm** (results in voxel size of 2.0 mm)

**Slice Thickness: 3.27 mm**

**Smoothing Filter:** NONE or 0.0 (for all filter options: loop filter, post-filter and z-axis filter)

All corrections 'On'

**Note:** Reconstruction offsets should be used to assure the head is fully in the in-plane FoV.

# PET Technical Procedures Manual

GE

Discovery STE - PET/CT scanners

(UCSD, UCSF)

## Acquisition Parameters:

### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

### Randoms Correction:

Singles: (not real-time subtraction)

## Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Reconstruction Method: 3D Iterative** (fully 3D Iter; not 3D FORE IR):  
**4 iterations; 20 subsets**

**Grid: 128 x 128**

**FOV: 256 mm** (results in voxel size of 2.0 mm)

**Slice Thickness: 3.27 mm**

**Smoothing Filter:** NONE or 0.0 (for all filter options: loop filter, post-filter and z-axis filter)

All corrections 'On'

**Note:** Reconstruction offsets should be used to assure the head is fully in the in-plane FoV.

# PET Technical Procedures Manual

GE

Discovery ST - PET/CT scanners

(UCI)

## Acquisition Parameters:

### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

### Randoms Correction:

Singles: (not real-time subtraction)

## Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Reconstruction Method: Iterative if available** (fully 3D Iter; not 3D FORE IR)

Only if fully iterative is not available, as in some older systems, is it ok to use 3D FORE IR.

**4 iterations; 21 subsets\*** (\* for older software versions having only 3D FORE IR, **24 subsets**)

**Grid: 128 x 128 x 47**

**FOV: 256 mm** (results in voxel size of 2.0 mm)

**Slice Thickness: 3.27 mm**

**Smoothing Filter:** NONE or 0.0 (for all filter options: loop filter, post-filter and z-axis filter)

All corrections 'On'

**Note:** Reconstruction offsets should be used to assure the head is fully in the in-plane FoV.



GE

## Omni Legend PET/CT scanner

# PET Technical Procedures Manual

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

**Randoms Correction:** Singles

### Reconstruction Parameters, phantom and all radiotracers:

**Reconstruction Method: Iterative:** VPHD (or VPFX): whichever you site normally uses, but NOT VPHDS or VPFXS

**8 iterations; 12 subsets**

**Grid: 256 x 256 matrix**

**FOV: 25.6 cm (results in 1.0 mm voxels)** Use Y-offset as needed to fit into this tighter FoV

**Smoothing: None** (there are 3 filters: in-plane, axial, and loop; all need to be **off**)

All corrections '**ON**'

# PET Technical Procedures Manual

## GE Signa - PET/MR scanners (Alabama, Stanford)

### Acquisition Parameters:

**Randoms Correction:** Singles

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Reconstruction Method: VPHD**  
**4 iterations; 16 subsets**

**Grid: 192 x 192 x 89**

**FOV: 256 mm** (results in voxel size of 1.333 mm)

**Slice Thickness: 2.79 mm**

**Smoothing Filter:** NONE or 0.0 (for all filter options: loop filter, post-filter and z-axis filter)

All corrections 'On'

**Note:** Reconstruction offsets should be used to assure the head is fully in the in-plane FoV.

## Philips Vereos - PET/CT scanners

# PET Technical Procedures Manual

(Vanderbilt)

## Acquisition Parameters:

### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

## Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Reconstruction Method: OSEM  
4 iterations; 21 subsets**

**Grid: 256 x 256 x 164**

**FOV: 256 mm** (results in voxel size of 1.0 mm) Slice Thickness: **1.0 mm**

**Smoothing:** Post-Gaussian filter should be 0 or OFF (no post-filtering)

The attenuation field should indicate “CTAC-SG” and the scatter field should indicate “SS-Simul”.

# PET Technical Procedures Manual

## Philips Ingenuity TF - PET/CT scanners (Penn)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

#### Acquisition Protocol: Brain

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

#### Reconstruction Method: Blob-OS-TF

**Grid: 128 x 128 x 90**

**FOV: 256 mm** (results in voxel size of 2.0 mm) Slice Thickness: **2.0 mm**

**Smoothing:** Set SMOOTHING parameter to '**SHARP**'

All other parameters should be set to defaults for the "**Brain**" protocol. All corrections '**On**'

For Blob-OS-TF reconstruction: The attenuation field should indicate "CTAC-SG" and the scatter field should indicate "SS-Simul".

# PET Technical Procedures Manual

## Philips Gemini TF - PET/CT scanners (Penn, Rush, Florida:Univ of Miami, Cleveland)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

#### Acquisition Protocol: Brain Protocol

### Reconstruction Parameters:

#### For reconstruction framing, see [Table 1](#):

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

#### Reconstruction Method: Iterative: LOR-Ramla

**Grid: 128 x 128 x 90**

**FOV: 256 mm** (results in voxel size of 2.0 mm)

**Slice Thickness: 2.0 mm**

**Smoothing:** Set SMOOTHING parameter to '**SHARP**'

All other parameters should be set to defaults for the "**Brain**" protocol. All corrections '**On**'

For LOR-Ramla reconstruction: The attenuation field should indicate "CTAC-SG" and the scatter field should indicate "SS-Simul".

# PET Technical Procedures Manual

## Siemens BioGraph Horizon - mCT TrueV PET/CT scanners (Wisconsin)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Method: Iterative: OSEM-3D (do NOT use TrueX reconstruction)**  
**4 iterations; 20 subsets**

**Grid: 360 x 360 x 81 (or 109 for the extended FoV model)**

**Zoom: 2.0** (results in voxel size of ~1.018 mm)

**Smoothing Filter: NONE** (All-pass or '0.0')

**Match CT slices: 'Off' or 'No'** (results in PET slice thickness of ~2.027 mm)

**All corrections 'On'**

I've seen 1 of 4 Horizon's that shows a slice thickness of 2.025 instead of 2.027 ??

# PET Technical Procedures Manual

**Siemens BioGraph mMR**  
(NYU, NY:Mount Sinai, Washington University)

## Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Method: Iterative: OSEM-3D (do NOT use TrueX reconstruction)**  
**4 iterations; 21 subsets**

**Grid: 344 x 344 x 127**

**Zoom: 2.0** (results in voxel size of 1.04313 mm)

**Smoothing Filter: NONE** (All-pass or '0.0')

**Match CT slices: 'Off' or 'No'** (results in PET slice thickness of 2.03125 mm)

**All corrections 'On'**

**Note for BioGraphs:** 1080, TruePoint (1093,1094), mCT, Horizon (possible mMR) To reconstruction into a matrix >256 your site has to have purchased the high-resolution option. Some have not, since it is actually pretty expensive. 256 is the best you can do. Resolution tends to be close to 0.5mm lower.

As you probably know, Siemens recon interpolated the sinogram projections down to the recon grid first, which losses MUCH resolution if you go down to 128. I know initially it was to speed up reconstruction, and then my (strong) guess is that at some point someone realized, when lets charge more for the ability to recon into large matrices and get the full resolution.

# PET Technical Procedures Manual

## Siemens BioGraph Vision PET/CT scanners (Indiana, Northwestern, New Mexico, UCSF, Washington University, Mayo)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Method: Iterative: OSEM-3D (do NOT use TrueX reconstruction)**  
**8 iterations; 5 subsets with ToF (time-of-flight)**  
**ToF must be used**

**Grid: 440 x 440 x 159**

**Zoom: 2.0** (results in voxel size of 0.825 mm)

**Smoothing Filter: NONE** (All-pass or '0.0')

**Match CT slices: 'Off' or 'No'** (results in PET slice thickness of ~1.645 mm)

**All corrections 'On'**



# PET Technical Procedures Manual

## Siemens BioGraph mCT - mCT TrueV PET/CT scanners

(UCDavis, New Mexico, Pittsburgh, Cleveland Clinic, USC, Florida:Univ of Florida, Columbia, Univ of Michigan, Yale)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

#### Scans and scan duration:

**LIST-MODE:** If your scanner has list-mode capability:

**NO LIST-MODE:** **ONLY** If your scanner does not have list-mode capability:

\*\*\* Note that reduce motion artifacts, two separate emission scans will be acquired as closely together as possible. Do not repeat CT scan.

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Method: Iterative: OSEM-3D (do NOT use TrueX or ToF reconstruction)**  
**4 iterations; 24 subsets**

**Grid: 400 x 400 x 81 (or 109 for TrueV model)**

**Zoom: 2.0** (results in voxel size of ~1.018 mm)

**Smoothing Filter: NONE** (All-pass or '0.0')

**Match CT **slices**: 'Off' or 'No'** (results in PET slice thickness of ~2.027 mm)

**All corrections 'On'**

# PET Technical Procedures Manual

## Siemens BioGraph HiRes – 81 slice PET/CT (Model 1080) (Florida: Mt Sinai)

### Acquisition Parameters:

#### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

#### PET acquisition:

**LIST-MODE:** If your scanner has list-mode capability

**NO LIST-MODE:** **ONLY** If your scanner does not have list-mode capability

\*\*\* Note that reduce motion artifacts, two separate emission scans will be acquired as closely together as possible. Do not repeat CT scan.

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Method: Iterative: FORE / OSEM-2D**  
**4 iterations; 16 subsets**

**Grid: 336 x 336 x 81**

**Zoom: 2.0** (results in voxel size of ~1.01567 mm)

**Smoothing Filter: NONE** (All-pass or '0.0')

**Match CT **slices**:** 'Off' or 'No' (results in PET slice thickness of ~2.00 mm)

**All corrections 'On'**

# PET Technical Procedures Manual

Siemens BioGraph TruePoint – and TruePoint TrueV  
PET/CT scanners (Models 1093, 1094)  
(Univ or Michigan, Washington University)

## Acquisition Parameters:

### CT scan:

Low mAs scan acquired shortly before emission. Leave enough time to start emission acquisition promptly at start time given immediately above. Scans should never be started early

### PET acquisition:

**LIST-MODE:** If your scanner has list-mode capability

**NO LIST-MODE:** **ONLY** If your scanner does not have list-mode capability

\*\*\* Note that reduce motion artifacts, two separate emission scans will be acquired as closely together as possible. Do not repeat CT scan.

## Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Method: Iterative: OSEM-3D (do NOT use TrueX reconstruction)**  
**4 iterations; 21 subsets**

**Grid: 336 x 336 x 81 (or 109 for TrueV (1094) model)**

**Zoom: 2.0** (results in voxel size of ~1.015 mm)

**Smoothing Filter: NONE** (All-pass or '0.0')

**Match CT **slices**: 'Off' or 'No'** (results in PET slice thickness of ~2.027 mm)

**All corrections 'On'**

# PET Technical Procedures Manual

## Siemens ECAT Exact HR+ (BGO) 63-slice scanners (Pittsburgh for PIB only)

### Acquisition Parameters:

#### Transmission scan:

**Five or six min 2-D scan** acquired immediately **post**-emission scan; process with **segmentation**.

### Reconstruction Parameters:

**For reconstruction framing, see [Table 1](#):**

Framing within the start-stop time in **Table 1** should be multiple 5 min frames (ex, for 30-60min start-stop time, 6x5min frames).

If dynamic data was acquired, uploading the reconstruction of the full dynamic dataset is encouraged as long as the frames within the start-stop time are reconstructed as multiple 5 min frames.

**Method: Iterative:** (FORE / OSEM-2D) **4 iterations; 16 subsets**

**Grid: 128 x 128 x 63**

**Brain Mode: ON**

**Zoom: 2.0**

**Smoothing Filter: NONE** (software version 7.2 says 'All Pass (Ramp)')

**Axial filtering: NONE** (software version 7.2 says 'Off')

All corrections '**On**'

# **Appendix B:**

## **IDA Uploader Introduction and User Registration (SCAN only)**

# PET Technical Procedures Manual

## INTRODUCTION

This document provides instructions for account registration and uploading images for the SCAN project in the Laboratory of Neuro Imaging's Image & Data Archive (IDA). For sites needing to upload site qualification scans, instructions for uploading to the SCANQUAL project are also included. The IDA utilizes a data de-identification process and encrypted file transmission to help ensure compliance with subject-privacy regulations.

## SYSTEM REQUIREMENTS

The IDA system requires the following:

- a computer with Internet access
- newer web browser software (IE/Edge, Firefox, Chrome, Safari)
- a valid user account with upload access for SCAN and/or SCANQUAL
- installation of the IDA Uploader application

## USER REGISTRATION *(Skip this step and go to [Appendix C](#) if you already have an IDA account)*

1. To register for a user account, go to the Image & Data Archive Log-In website (<https://ida.loni.usc.edu>) and select “Instant free signup” in the top right corner.



2. Complete New account registration (3 steps).
  - a. Enter your email address and select “CONTINUE”. A security code will be emailed to you.

# PET Technical Procedures Manual

## New account registration - Step 2 of 3

To complete the email verification process, enter the security code sent to your email address below.

- b. Enter the security code sent to your email and click “CONTINUE”.

- c. Complete the New account registration form and click “REGISTER”. A link to set your password will be emailed to you with subject line “Welcome to the LONI Image & Data Archive”.

## New account registration - Step 3 of 3

Click REGISTER to complete the new account registration process. You will receive an email containing a link to the page to set your password.

- d. Create a password and click “CONTINUE”.

## Create Password

3.

Your account is created, if you need upload access to SCAN or SCANQUAL, please send an email to the appropriate contact following the steps below:

For access to SCAN or SCANQUAL, email [data.coordinator@loni.usc.edu](mailto:data.coordinator@loni.usc.edu)

- a. Enter “SCAN Upload Access Request” or “SCANQUAL Upload Access Request” in the subject line of your email.
- b. Provide the email address you used when creating your account, your site name and site number in your email request.

## **PET Technical Procedures Manual**

- c. You will receive an email when your account access has been set, generally within one working day.

# **Appendix C: Obtaining/Installing IDA-uploader (SCAN only)**



# PET Technical Procedures Manual

## OBTAINING AND INSTALLING THE IDA-UPLOADER

- a. [Installing the IDA-Uploader for Windows](#)
- b. [Installing the IDA-Uploader for Mac](#)
- c. [Installing the IDA-Uploader for Linux](#)

### a. *Installing the IDA-Uploader for Windows*

1. Log in to the IDA and select SCAN from the PROJECTS menu.
2. Select the ARCHIVE Menu option.
3. Select your operating system (Windows 32-bit or 64-bit) from the dropdown menu.
4. Click “Download”.

HOME	SCAN @LONI	DOWNLOAD	SEARCH	ARCHIVE	MANAGE	PROJECTS	SUPPORT
------	------------	----------	--------	---------	--------	----------	---------

The upload process involves two basic steps:

- 1 De-identify file metadata by replacing any fields that identify the subject, such as Patient Name and ID.
- 2 Transmit files securely from the local site to LONI.

### IDA Uploader Application



You will need to launch the Uploader application from your computer to upload SCAN data. Choose your operating system and download the application below.

Operating System


Windows 64-bit ▾

Download


5. Open the application by clicking on the download in your browser or by locating the IDA-Uploader-2.0.msi application in the Downloads section of your File Explorer.

# PET Technical Procedures Manual

6. You will be taken to the IDA-Uploader Setup Wizard – click “Next” to continue.



7. Choose your destination folder and click “Next”.




8. Click “Install”.


# PET Technical Procedures Manual

Note: If you

receive a popup asking for permission to run the application, click “yes”.



9. Click “Finish” to exit the Setup Wizard.



10. Locate the application on your device by typing in “IDA-Uploader” in the Windows Start Menu.

## ***b. Installing IDA-Uploader for Mac***

1. Log in to the IDA and select SCAN from the PROJECTS menu.

# PET Technical Procedures Manual

2. Click ARCHIVE from the Menu.
3. Select your operating system (Mac) from the dropdown menu.
4. Click “Download”.

HOME	SCAN @LONI	DOWNLOAD	SEARCH	ARCHIVE	MANAGE	PROJECTS	SUPPORT
------	------------	----------	--------	---------	--------	----------	---------

The upload process involves two basic steps:

- 1 De-identify file metadata by replacing any fields that identify the subject, such as Patient Name and ID.
- 2 Transmit files securely from the local site to LONI.

## IDA Uploader Application



You will need to launch the Uploader application from your computer to upload SCAN data. Choose your operating system and download the application below.


Operating System

Mac

Download


5. Open your Finder application, click Downloads, right-click ‘IDA-Uploader-2.0.pkg’ and select “Open”.

NOTE: If you try to install by double-clicking on ‘IDA-Uploader-2.0.pkg’, you may receive a message that it cannot be opened. Please ensure that you follow the instructions above to successfully open.




# PET Technical Procedures Manual

6. An Install IDA-Uploader window will appear. Click “Continue”.




7. Next, choose a destination for the installation. Click “Continue”.




# PET Technical Procedures Manual

8. For Installation Type, you can review the details of the installation. Click “Install”.




9. During the installation, you will need to enter the User Name and Password for the computer. Click “Install Software”.



# PET Technical Procedures Manual

10. Once the installation is complete, the window will provide a summary of a successful installation. Click “Close”.

NOTE: A window will appear to ask “Do you want to move the “IDA-Uploader” installer to the trash?”. You can choose Keep or Move to Trash.



11. Locate the application in your computer by opening your Finder application. Click on Applications and double-click IDA-Uploader.

# PET Technical Procedures Manual

## c. Installing IDA-Uploader for Linux

1. Log in to the IDA and select SCAN from the PROJECTS menu.
2. Select the ARCHIVE Menu option.
3. Select your operating system (Linux) from the dropdown menu.
4. Click “Download”.



The upload process involves two basic steps:

1

De-identify file metadata by replacing any fields that identify the subject, such as Patient Name and ID.

2

Transmit files securely from the local site to LONI.

### IDA Uploader Application



You will need to launch the Uploader application from your computer to upload SCAN data. Choose your operating system and download the application below.

Operating System

Linux

Download

5. Please visit <https://www.oracle.com/java/technologies/javase-downloads.html> to download the latest Oracle JDK. NOTE: A minimum version of 15.0.1 is required to run IDA-Uploader-2.0.jar
6. Choose DEB, PRM or the compressed archive depending on your Linux Distribution.
7. Once you have completed the installation, open the terminal and run: `java -jar IDA-Uploader-2.0.jar`



## **Appendix D:** How to Upload PET Data to SCAN using IDA-uploader (SCAN only)


# PET Technical Procedures Manual

## UPLOADING (ARCHIVING) TO SCAN

1. Open the IDA-Uploader application.
2. Enter your email and password, then click “Log In”.




3. Select SCAN from the Project dropdown menu.



# PET Technical Procedures Manual

4. Then, select your site from the Site dropdown menu. Click “Continue”.



5. Enter the participants identifiers into the Subject ID field using the convention PTID+NACC ID, where the PTID and NACC ID are separated by the '+' character. As a rule, the user should enter the complete PTID and NACC ID and follow the formatting convention used for UDS data submission when applicable. Either upper or lower case may be used for any alphabetic characters. If the given PTID and NACC ID pair match the record present in the NACC database, the upload will be accepted. If not, the user will receive an error message indicating that the PTID and NACC ID combination are not valid, and the upload will not be allowed to continue.

**PLEASE NOTE: To upload data you must enter both the PTID and the NACC ID. To access your center's list of PTID and NACC ID pairs, you may work with your center UDS data manager to utilize the "PTID to NACC ID Map" tool available via the NACC portal:**  
<https://www.alz.washington.edu/MEMBER/portal>.

Participant identifier terminology:

- PTID (may be referred to as Subject ID or local ID). This is the ADRC-managed participant ID which is a sequence of up to 10 characters. Formats vary by each center; this is the ID used for UDS (uniform dataset) identification internally.
- NACC ID: This is the NACC-managed participant ID. It is a string with the prefix 'NACC' followed by 6 digits.

6. Click “Browse” to select the Source Directory. Then click “Upload”.

NOTE: The Source Directory is the directory containing the files to be uploaded. If your Source Directory contains subdirectories, choose to include/exclude those files by checking “Search subdirectories”.

# PET Technical Procedures Manual

IDA Uploader 2.0

Image & Data Archive  
Powered by LONI

**IDA Uploader**  
Version 2.0

Reference ID:  
SCAN-ACD10-1612557673311

Select Files De-Identify and Upload Complete Upload in Browser

Enter required information below and click Upload to begin the process.

**Site**  
Mayo Clinic (ADC10)

**Subject ID**  
De-identifier to replace patient ID. Maximum of 10 characters allowed.  
Enter Subject ID

**Source Directory**  
Location of files to de-identify and upload.  
Click Browse to select directory Browse


☐ Search subdirectories

Need help? Email [ida@loni.usc.edu](mailto:ida@loni.usc.edu).

Cancel Upload


7. You will be able to see the progress of your upload in the De-identify and Upload section.

# PET Technical Procedures Manual



8. Once the files are de-identified and transferred to the IDA, you will need to complete the upload in your web browser.

# PET Technical Procedures Manual



9. Your web browser should automatically open a new "Log In to Continue" page. Enter your IDA email and password. Click "Log In".

# PET Technical Procedures Manual




Image & Data Archive  
Powered by LONI

Log In to Continue

Email Address

Password


[Forgot password? Click here](#)

[Log In](#) [Sign Up](#)

© 2003 - 2023 LONI. All rights reserved.

NOTE: If your web browser does not automatically open the “Log In to Continue” page, please click “Complete Upload” and you can copy-paste the link in your browser.

# PET Technical Procedures Manual




10. In the “Your Upload is Not Complete” page you can review the details of your upload.



# PET Technical Procedures Manual


[HOME](#) [SCAN @LONI](#) [DOWNLOAD](#) [SEARCH](#) [ARCHIVE](#) [MANAGE](#) [PROJECTS](#) [SUPPORT](#)

 **Your Upload is Not Complete for TEST0123**


Additional information is required to complete the archiving process. Please review and complete the details below. All fields are required.

Uploader Email: @loni.usc.edu  
Site: ADC10  
Upload Date: February 11, 2021 9:43 AM PST

Subject Details: Participant

 MRI Images (1)

Description	Series Date	Image	Delete
MPRAGE GRAPPA	November 16, 2012 10:24 AM	<a href="#">View</a>	<a href="#">X</a>

 PET Images (1)

Description	Series Date	Image	Metadata	Delete
Downs (128x128,3mm)	May 22, 2017 3:32 PM	<a href="#">View</a>	<a href="#">Complete Form</a>	<a href="#">X</a>

After completing click Finish Upload.

[X Remove All](#) [✓ Finish Upload](#)

- a. For PET image uploads, additional information is required in order to complete the upload.
  1. In the Metadata section, click “Complete Form”.

NOTE: Fields marked (\*Required) are required to complete the upload.
  2. Once the information is entered in the form, click “Update”.

# PET Technical Procedures Manual





## Complete Metadata Form

Date July 28, 2022 5:51 PM

Visit

Description Florbetaben.4x5min.3D.440x.8i5sTOF.2z.AllPass


Scan Type	
<b>Scan Date</b> * Required	 September 23, 2022
<b>Tracer</b> * Required	Select One 
<b>PET tracer dose assay (mCi)</b> * Required	PET tracer dose assay (mC)
<b>Time of tracer dose assay (24-hour clock)</b> * Required	HH : MM
<b>Time of injection (24-hour clock)</b> * Required	HH : MM
<b>Emission start time (24-hour clock)</b> * Required	HH : MM
Data Transfer	
<b>Additional comments</b>	<div>Please explain in 384 characters or less</div> <div>384</div>

Cancel

Update

# PET Technical Procedures Manual

NOTE: If any data was mistakenly uploaded, you can click the “X” in the Delete column and it will be deleted from your upload. A window will appear to confirm the deletion. Click “Delete” to confirm file deletion or “Cancel”.





# PET Technical Procedures Manual

13. Once the upload is processed, a summary page for the upload will be displayed. Please note that you have the option of clicking “Download CSV Files” to keep a record of the upload.

Summary for TEST0123 SCAN\_ADC10-1613066168247

A summary of your upload is below. You can copy the reference URL or download this summary as a CSV file for future reference.

**Download CSV Files**

Email address of uploader	Site	Upload Date	Archive Date	Research Group	Reference URL
@loni.usc.edu	ADC10	February 11, 2021 10:42 AM PST	February 11, 2021 9:57 AM PST	Participant	Copy URL

**MRI Images**

Description	Series Date	Image
MPRAGE GRAPPA	November 16, 2012 10:24 AM	<a href="#">View</a>

**PET Images**

Description	Series Date	Image	Metadata
Downs (128x128,3mm)	May 22, 2017 3:32 PM	<a href="#">View</a>	<a href="#">View</a>

14. You can close the IDA-Uploader application or to upload images for another subject, click the “Upload More” button in the IDA-Uploader application.

# **Appendix E:**

## **How to Upload PET Hoffman Phantom Data to SCANQUAL (SCAN only)**

# PET Technical Procedures Manual


## UPLOADING (ARCHIVING) TO SCANQUAL

1. Open the IDA-Uploader application.
2. Enter your email and password, then click “Log In”.



The screenshot shows the IDA Uploader 2.0 application window. The title bar reads "IDA Uploader 2.0". The main content area features the "Image & Data Archive Powered by LONI" logo at the top, followed by the text "Login to the Image and Data Archive". Below this, there is a login form with fields for "Email" and "Password". The "Email" field is a simple text input. The "Password" field is a text input with a placeholder "Enter Password". To the left of the password field is a circular icon containing a padlock. Below the password field, there is a link "Forgot Password? Click here" and a "Log In" button, which is highlighted with a red rectangle. To the right of the "Log In" button is a "Sign Up" button.

3. Select SCANQUAL from the dropdown menu. Once you select SCANQUAL, a dropdown with the sites will appear. Select your site and click Continue.



The screenshot shows the IDA Uploader 2.0 application window. The title bar reads "IDA Uploader 2.0". The main content area features the "Image & Data Archive Powered by LONI" logo at the top, followed by the text "Choose a Project/Site". Below this, there is a form with a "Project" label and a dropdown menu. The dropdown menu is highlighted with a red rectangle and contains the text "Select Project". Below the dropdown menu is a "Continue" button.

# PET Technical Procedures Manual

IDA Uploader

Single Archive

Batch Archive

4. Click "Single Archive" or "Batch Archive".

## a. Single Archive

Use the Single Archive process to upload one or more files from a single subject.

1. After clicking "Single Archive" the De-Identification page will be displayed.

Project SCANQUAL@ADC33 ☐ Bypass validation steps

Select Data Type ☒ Original ☐ XML

Subject ID:  Max. 10 characters allowed  
Identifier to replace Patient ID

Source Directory:    
Location of original files

Target Directory:    
Location for target files

NOTE: The Source Directory may contain multiple data directories and data files for the same subject and visit.

☐ Record diagnostics to file

- a. Click the type of data being uploaded – in this case, "Original".

- b. Enter the Phantom/Volunteer ID in the Subject ID field.

**Phantom:** 1\_P\_0001 (siteID\_P\_Number)

**Volunteer:** 1\_V\_0001 (siteID\_V\_Number)

Assign the Phantom/Volunteer ID in the Subject ID field as: "siteID\_P\_xxxx" or "siteID\_V\_xxxx", where siteID is your 1-3 digit Site ID and xxxx should count up for your site starting with 0001. Do not add leading zeros for siteID.



For

# PET Technical Procedures Manual

example, if your Site ID is ADC 13, you should fill in “13\_P\_0002” for the *second* phantom scan to be uploaded.

- c. Select the **Source Directory** in which the original files are located.
- d. Select the **Target Directory** for de-identified files to be written to.
- e. Click “CONTINUE” to begin the de-identification process.

## 2. On the Verify and Submit page

- a. Deselect any image you do not want to be archived (if any) by unchecking the Selected checkbox.
- b. Click “SUBMIT” to begin the transmission process.

Note: This is not a feature during Batch Archive. Once the transmission has begun, a progress bar will show the status of the upload.

The screenshot shows the 'IDA Uploader' window. At the top, there are fields for 'Subject ID' (129\_V\_9898), 'Sequence Name' (MPRAGE\_3dtferrepeat), 'Number of Images' (1), and 'Selected' (1). Below these fields are two buttons: 'DISCARD' and 'SUBMIT'. The 'SUBMIT' button is highlighted with a red rectangular box. To the right of the buttons is a checkbox labeled 'Compress files before transmitting', which is checked. Below the buttons is a section titled 'REVIEW DE-IDENTIFIED HEADER INFORMATION:'. This section contains a list of DICOM file attributes and their values, including 'Series Description', 'Series ID', 'Metadata for DICOM file', and a table of tags and values.

Tag	Tag Description	Tag Value
00020000	Group Length	180
00020001	File Meta Information Version	<BYTE>
00020002	Media Storage SOP Class UID	1.2.840.10008.5.1.4.1.1.4
00020003	Media Storage SOP Instance UID	2.16.124.113543.6006.99.07932364392440536066
00020010	Transfer Syntax UID	1.2.840.10008.1.2.2
00020012	Implementation Class UID	2.16.124.113543.6006.99.7256096479968646091
00080000	Group Length	1064
00080005	Specific Character Set	ISO IR 100
00080008	Image Type	ORIGINAL PRIMARY M_FFE M_FFE
00080013	Instance Creation Time	120651
00080014	Instance Creator UID	2.16.124.113543.6006.99.08781057918761804044
00080016	SOP Class UID	1.2.840.10008.5.1.4.1.1.4
00080018	SOP Instance UID	2.16.124.113543.6006.99.07932364392440536066
00080020	Study Date	20080603
00080021	Series Date	20080603
00080022	Acquisition Date	20080603
00080023	Image Date	20080603
00080030	Study Time	112614
00080031	Series Time	114448.98000
00080032	Acquisition Time	114448.98
00080033	Image Time	114448.98
00080060	Modality	MR
00080070	Manufacturer	Philips Medical Systems
00080080	Institution Name	Dartmouth College
00080090	Referring Physician's Name	
00081030	Study Description	ADNI
0008103E	Series Description	MPRAGE_3dtferrepeat

### b. Batch Archive

The Batch Archive process is similar to Single Archive, except that multiple subjects and image series can be submitted in a batch. Batches can be of the same or different modalities. However, users cannot review the results of the de-identification process prior to the batch upload.

1. Proceed to follow the De-identification steps in the Single Archive section.
2. The Batch Archive will skip the Verify and Submit step that is available in Single Archive, and direct you to the Image Database Batch Queue page.
3. Click “ADD MORE” to add more images to the Batch. Repeat this process until you have added everything you intend to archive.
4. Click “SUBMIT” to begin both the de-identification and transmission processes.

# PET Technical Procedures Manual

IDA Uploader							
Subject	Data Type	Research Group	Source		Status	Date	Remove
131_P_5555	Original	Phantom	/Users/t	/Desktop...	Queued	2/08/21	<a href="#">remove</a>
131_P_9090	Original	Phantom	/Users/t	/Desktop...	Queued	2/08/21	<a href="#">remove</a>
131_V_3333	Original	Volunteer	/Users/t	/Desktop...	Queued	2/08/21	<a href="#">remove</a>
131_V_3333	Original	Volunteer	/Users/t	/Desktop...	Queued	2/08/21	<a href="#">remove</a>
131_V_4444	Original	Volunteer	/Users/t	/Desktop...	Queued	2/08/21	<a href="#">remove</a>
131_P_6767	Original	Phantom	/Users/t	/Desktop...	Archived	2/08/21	<a href="#">remove</a>

Text

---

## **Appendix F:** **SCAN PET Upload Form** **(SCAN only)**

Site: \_\_\_\_\_

PTID or Subject ID : \_\_\_\_\_

Scan date (mm/dd/yyyy): \_\_\_\_\_

PET tracer: \_\_\_\_\_

PET tracer dose assay (mCi): \_\_\_\_\_

Time of tracer dose assay (military time hh:mm): \_\_\_\_\_

Time of tracer injection (military time hh:mm): \_\_\_\_\_

Emission start time (military time hh:mm): \_\_\_\_\_

#### Possible PET tracers:

PIB

Florbetapir

Florbetaben

NAV4694

Flortaucipir

MK6240

PI2620

GTP1

FDG